



X-RAY MACHINE 200MA

Control

- * 100 KVP 200 mA
- * Normal frequency 220 V 50 HZ
- * Contains the keys to control the following factors:
 - * KV: From 40 to 100 kVP. Control up to 1 kVP
 - * mA: from 50 to 200 mA
 - * Time: From 0.01 to 3 seconds

X-Ray Tube

JX51-20 ,40/125 made in China 200mA

Tube Stand

The tube holder moves vertically to change the shooting distance and moves horizontally on the longitudinal axis of the table

4 Way Floating Table

- * Imaging table with mechanical compass
- * Moving table in 4 directions and weights up to 180 kg.
- * Sizes of the cassette place from 8 x 10 inches up to 17 x 17 inches
- * Fine grid x-ray filter orthogonal lines 103 lines ratio 1:10
- * Installation of the imaging table does not require any installation ceiling or walls or any modification in the room

Chest Stand

- Contains a cassette cabinet that moves up and down from ground level up to 200 cm in sliding motion system
- * And moving weights, with mounting brakes (manual).
- * Variable cassette holder sizes from 8 x 10 inches up to 17 x 17 inches
- * Contains linear ray filter orthogonal lines 103 line ratio 1:10 made in Korea



