



# HKN-9010

## INFANT RADIANT WARMER

### Brief Introduction

- Three control Modes: Pre-warm Mode, Manual Mode and Baby Mode control by micro-computer;
- Set Temperature and Baby Temperature can be displayed separately;
- Horizontal angle of warmer module and inclination of the bassinet can be adjusted;
- The panels surrounded bassinet can be turned outward or disassembled;
- X-ray cassette under the bassinet;
- APGAR timer;
- Multiple failure alarm indications;
- RS232 connector for data output;
- Disposable skin temperature sensor and SpO<sub>2</sub> are optional.

### Standard Configuration

Main body (including the warmer module, bassinet, main column, mobile stand and temperature controller), I.V. Pole and tray.

### Specification

Power requirement:	AC220V~230V/50Hz or AC110~120V/50~60Hz, 600VA
Control Mode:	Pre-warm Mode, Manual Mode and Baby Mode are controlled by micro-computer
Baby Mode temperature control range:	34.5~37.5℃
Temperature sensor display range:	5℃~65℃
Accuracy of skin temperature sensor:	±0.2℃
Temperature uniformity of mattress:	≤2℃
Angle of heater head:	0°~90° any angel for two ways
Inclination of the bassinet:	three tilt position
APGAR timer:	audible and visible alarm when the device runs to 50" ~ 1' , 4' 50" ~ 5' , 9' 50" ~ 10'
Failure alarm:	over temperature alarm, deviation alarm, sensor failure alarm, power failure alarm, setting alarm, checking alarm and so on.

### Environmental Requirement

Operating range: 18℃~30℃  
Ambient air velocity rate: <0.3m/s

### Product Packaging

Package: The whole unit is packed into one carton  
Dimension: 1220mm x 670mm x 1100mm  
Gross Weight: 100kg



Stock code [300314]



Ningbo David Medical Device Co.,Ltd

Sales center: #100 Jinghua Road, High-tech Industrial Development Zone, Ningbo, China  
Tel: 0086-574-8780 0008 8780 0007 /Fax: 0086-574-8780 3714  
E-mail: sales@nbdavid.com /Http: //www.chinadavid.cn /P.C.: 315040