



Model	JH20-1B
Checking methods	Compare the blue and green light
Indication	Three-figure LED digital display(one decimal digits)
Readable	Test the bilirubin value directly by the forehead of baby
Average counting	No
Test deviation	0~15±1 16~25±1.5
Light resource	Xenon-flashing light
Power supply	Nickel xenon battery group 4.8V
Prepare time of start-up	<10s
Charger	Input 220V 50Hz Output 6V 0.3ADC
Checkout	White screen indicates 00 or 01 Yellow screen indicates 20.0±1
weight	245g

