

#### Control

- \* 110 KVP 500 mA
- \* Normal frequency 220 V 50 HZ
- \* Contains the keys to control the following factors:
  - \* KV: From 40 to 110 kVP. Control up to 1 kVP
  - \* mA: from 50 to 500 mA
  - \* Time: From 0.01 to 5 seconds

## X-Ray Tube

Canon model E7239X made in Japan 500 mA

# **Tube Stand**

The tube holder moves vertically to change the shooting distance and moves horizontally on the longitudinal axis of the table

# 4 Way Floating Table

- \* Imaging table with mechanical compass
- \* Moving table in 4 directions and weights up to 180 kg.
- \* Sizes of the cassette place from 8 x 10 inches up to 17 x 17 inches
- \* Fine grid x-ray filter orthogonal lines 103 lines ratio 1:10
- \* Installation of the imaging table does not require any installation ceiling or walls or any modification in the room

## **Chest Stand**

- \* Contains a cassette cabinet that moves up and down from ground level up to 200 cm in sliding motion system
- \* And moving weights, with mounting brakes (manual).
- \* Variable cassette holder sizes from 8 x 10 inches up to 17 x 17 inches
- \* Contains linear ray filter orthogonal lines 103 line ratio 1:10 made in Korea



